|                                                                            |           |
|----------------------------------------------------------------------------|-----------|
| <b>Instructions for Use with Pressure Thermoforming Equipment</b>          | <b>2</b>  |
| <b>Mode d'emploi pour les dispositifs formés par pression</b>              | <b>5</b>  |
| <b>Istruzioni per l'uso con macchine termoformatrici a pressione</b>       | <b>8</b>  |
| <b>Gebrauchsanweisung für die Verarbeitung mit Druck-Tiefziehgeräten</b>   | <b>11</b> |
| <b>Instrucciones de uso con equipos de termoconformado por presión</b>     | <b>14</b> |
| <b>Instruções de uso com equipamentos termoformadores de pressão</b>       | <b>17</b> |
| <b>Інструкція із застосування для обладнання тискового термоформування</b> | <b>20</b> |

## Zendura VIVA® Thermoformable Sheet Materials

Zendura VIVA single sheet and multipacks • REF# 9212, 9251, 9284, 9287, 9252

Zendura VIVA 20-sheet box • REF# 9212-20, 9284-20, 9287-20

### Instructions for Use with Pressure Thermoforming Equipment

#### Intended Use

Zendura Thermoformable Materials are intended for the fabrication of orthodontic and dental appliances designed and planned by a registered orthodontist or registered dental lab technician.

#### Indications for Use

Zendura Thermoformable Materials allow the fabrication of orthodontic and dental appliances such as aligners and retainers.

|                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <b>Contraindications</b><br>This material is contraindicated for patients with a history of allergic reaction to plastics.                                                                                                                                                                                                                                                                                                                                                                                           |
|    | <b>Warning:</b> Keep Zendura thermoformable material away from excessive heat or ignition source. Provide suitable exhaust ventilation when thermoforming or trimming, as it may potentially form combustible dust.                                                                                                                                                                                                                                                                                                  |
|    | <b>Warning:</b> Excessive heating time during thermoforming will cause plastic to sag with possible risk of injury to operator.                                                                                                                                                                                                                                                                                                                                                                                      |
|    | <b>Caution:</b> Zendura is provided pre-dried in moisture barrier packaging and should be thermoformed <b><u>within 15 minutes</u></b> of opening its protective foil pouch. This time can be reduced in a humid environment. Delayed thermoforming may cause bubbling and other defects that may render the final dental appliance unsuitable for use.                                                                                                                                                              |
|    | <b>Caution:</b> Failure to properly prepare the model, thermoform, trim, or clean the appliance, or improper handling in general, may result in defects, cracking, or breakage that may create sharp edges that could cause injury, loose pieces that could be aspirated, and/or render the dental appliance unsuitable for use.                                                                                                                                                                                     |
|  | <b>Caution:</b> Store unopened Zendura packages in a cool and dry place in their original packaging. Do not use Zendura thermoformable material if the moisture barrier bag sealing is damaged or compromised. Do not use Zendura thermoformable material after the expiration date indicated on the packaging. Product performance may be affected.<br>Dispose of material and used dental appliances in accordance with applicable local regulations. Used dental appliances may be considered biohazardous waste. |
|  | <b>Caution:</b> Zendura thermoformable material is intended for single use only. Reuse or repeated thermoforming may compromise the material's mechanical properties and device performance.                                                                                                                                                                                                                                                                                                                         |

**Important:** This product is for use by a professional user, according to the instructions for use.

**Important:** Users and/or patients should report serious incidents with medical device information to the manufacturer and to the national competent authority.

### Model Preparation

1. The dental model should be well cured and dry. If using 3D printed models, follow all post-processing steps recommended for the resin.
2. The model should be trimmed to be about 11 mm – 15 mm (0.4"-0.6") high. For best results the model should be shaped in the shape of an arch without the palate or tongue space filled-in.
3. Fill in any holes or undercuts on the model in the model design software or with block out material and let dry or cure.
4. If needed, coat the model with model release agent (foodgrade silicone or non-stick cooking spray) to prevent the appliance from sticking to the model and let dry.

### Thermoforming

1. Turn on the thermoformer and enter correct code or adjust timer for Zendura sheet
2. (see chart with recommended configuration settings).
3. Thermoform using a pressure of 4 bar (60 psi) or above. A higher pressure is preferable.
4. Place the model in the center of the thermoforming platform.
5. Preheat the thermoformer's heating element before placing the model on the platform.
6. Open the foil pouch. Place the Zendura sheet in the thermoforming frame and secure it.
7. Swing the heating element over the Zendura sheet to begin heating until the sheet is ready to be formed
8. Heat material for the recommended time.
9. Move the heating element away from the sheet.
10. Form the Zendura sheet over the model with pressurized air and allow the part to cool before handling.
11. Remove the thermoformed model from the frame.
12. Clean according to instructions given further below

### Trimming Instructions

1. Trim and finish the appliance in a well-ventilated environment. It is recommended to wear a mask and trim with a vacuum unit to reduce dust. Do not use hot knife for finishing the appliance. Polish the edges.
2. Excess material can be removed using utility scissors to make the model easier to handle.
3. Cut Zendura appliance off the dental model using a serrated trimming wheel or laboratory bur. Operate the trimming wheel or laboratory bur at around 35,000 rpm (medium speed) and cut as close to the finished outline of the appliance as possible.
4. Note: It can be helpful to draw the finish line margin on the thermoformed material as a guide prior to trimming.
5. Remove appliance from the 3D model, trim off extra parts with crown and bridge scissors or coarse rotary polishing wheel.
6. Polish the edges of Zendura appliance with a rotary polishing wheel at a speed of about 10,000 rpm. Medium and fine polishing wheels are recommended.
7. Clean the appliance according to the instructions provided further below.

### Suggested Zendura VIVA Pressure thermoforming settings

- If the plastic does not form well to the model, add 5 seconds to heating time until conformation is acceptable.
- If the plastic forms webs or folds, reduce heating time until webbing does not occur.

### Suggested Thermoforming Settings

**Note:** Times are guidelines only, as heating times vary between different thermoformers.

|                 | Temperature     | Heating Time | Pressure | Code |
|-----------------|-----------------|--------------|----------|------|
| MiniSTAR® S     | 220 °C          | 35 sec       | ≤4 Bar   | 132  |
| Biostar®        | 220 °C          | 45 sec       | ≥4 Bar   | 152  |
| Druformat Scan  |                 | 1 min 15 sec | ≥4 Bar   |      |
| Erkodent® Units | 190 °C - 205 °C |              |          |      |

*MiniSTAR, MiniSTAR S and Biostar are trademarks of Scheu Dental Technology. Druformat Scan is trademark of Dreve Dentamid. Erkodent is trademark of ERKODENT Erich Kopp GmbH.*

**Cleaning**

Wash Zendura appliance with mild soap and rinse well with water. Do not use hot water, vinegar, organic solvents or alcohol.

**Storage and Disposal**

Store unopened and sealed Zendura thermoformable material in a cool and dry place in their original packaging. Disposal should be handled in an environmentally sustainable manner according to local regulations. Used dental appliances may be considered a biohazard. Follow your country-specific laws, directives, standards and guidelines for the disposal of used devices.

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\*U.S. Pat# 10,870,263 & 10,549,511; International patents pending. ZENDURA is a registered trademark of Bay Materials, LLC.
© Bay Materials, LLC. All rights reserved.
*MiniSTAR, MiniSTAR S and Biostar are trademarks of Scheu Dental Technology. Druformat Scan est une marque déposée de Dreve Dentamid.*
Emergo Asia Pacific Pty Ltd
Bay Materials, LLC
Etkon GmbH